

File D:\DATA\

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN

Gradient: 5%-95% B in 1.5min


Flow Rate: 2.0ml/min

Column: XBridge C18,4.6*50mm,3.5um


Oven Temperature: 40 C

MASS Range:110-800


*ELS1 A, Voltage




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ELS1 A, Voltage

RetTim	Width	Area	Height	Area%
1.80	0.01	416.46	477.11	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.80	0.01	416.46	477.11	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.80	0.01	828.44	946.85	100.00

Ret. Time: 1.80

<<<< POSITIVE SPECTRA >>>>

