

YX-118292-40-3
YX-118292-40-3

9.02
9.02

8.27
8.26
8.25
8.24

7.73
7.71
7.63
7.62

7.26
7.26
7.24
7.24

7.10
7.08

4.35
4.33
4.31
4.30

3.30
3.04
3.03
3.02
3.01

2.47
2.47
2.46

1.87
1.86
1.84

1.32
1.31
1.29
1.26

Parameter	Value
1 Spectrometer	vnmrs
2 Solvent	dmso
3 Temperature	25.0
4 Pulse Sequence	s2pul
5 Number of Scans	8
6 Receiver Gain	34
7 Relaxation Delay	1.0000
8 Pulse Width	0.0000
9 Acquisition Time	4.0630
10 Acquisition Date	2017-07-27T16:13:09
11 Spectrometer Frequency	399.74
12 Spectral Width	8012.8
13 Lowest Frequency	-1608.0
14 Nucleus	1H
15 Acquired Size	32556
16 Spectral Size	65536

