

265121-04-8 (HX201707073102)-HNMR-20170731  
YX-165121-04-8

| Parameter                 | Value               |
|---------------------------|---------------------|
| 1 Spectrometer            | vnmrs               |
| 2 Solvent                 | dmsd                |
| 3 Temperature             | 25.0                |
| 4 Pulse Sequence          | s2pul               |
| 5 Number of Scans         | 8                   |
| 6 Receiver Gain           | 36                  |
| 7 Relaxation Delay        | 1.0000              |
| 8 Pulse Width             | 0.0000              |
| 9 Acquisition Time        | 4.0630              |
| 10 Acquisition Date       | 2017-07-31T16:23:33 |
| 11 Spectrometer Frequency | 399.74              |
| 12 Spectral Width         | 8012.8              |
| 13 Lowest Frequency       | -1608.0             |
| 14 Nucleus                | 1H                  |
| 15 Acquired Size          | 32556               |
| 16 Spectral Size          | 65536               |

