



Parameter	Value
1 Solvent	dms
2 Temperature	80.0
3 Pulse Sequence	s2pul
4 Number of Scans	8
5 Receiver Gain	32
6 Relaxation Delay	1.0000
7 Pulse Width	0.0000
8 Acquisition Time	4.0630
9 Acquisition Date	2017-08-21T10:15:21
10 Modification Date	2017-08-22T01:16:27
11 Spectrometer Frequency	399.74
12 Spectral Width	8012.8
13 Lowest Frequency	-1608.0
14 Nucleus	1H
15 Acquired Size	32556
16 Spectral Size	65536