

Parameter	Value
1 Solvent	dmsd
2 Temperature	25.0
3 Pulse Sequence	s2pul
4 Number of Scans	8
5 Receiver Gain	44
6 Relaxation Delay	1.0000
7 Pulse Width	0.0000
8 Acquisition Time	4.0630
9 Acquisition Date	2017-12-21T14:26:36
10 Modification Date	2017-12-22T06:27:36
11 Spectrometer Frequency	399.74
12 Spectral Width	8012.8
13 Lowest Frequency	-1608.0
14 Nucleus	1H
15 Acquired Size	32556
16 Spectral Size	65536

