

Parameter	Value
1 Solvent	DMSO
2 Temperature	298.1
3 Pulse Sequence	zg30
4 Number of Scans	8
5 Receiver Gain	94
6 Relaxation Delay	1.0000
7 Pulse Width	8.0000
8 Acquisition Time	3.2768
9 Spectrometer Frequency	500.13
10 Spectral Width	10000.0
11 Lowest Frequency	-1913.2
12 Nucleus	<sup>1</sup> H
13 Acquired Size	32768
14 Spectral Size	65536

