

| Parameter                | Value          |
|--------------------------|----------------|
| 1 Solvent                | DMSO           |
| 2 Temperature            | 298.1          |
| 3 Pulse Sequence         | zg30           |
| 4 Number of Scans        | 8              |
| 5 Receiver Gain          | 167            |
| 6 Relaxation Delay       | 1.0000         |
| 7 Pulse Width            | 8.0000         |
| 8 Acquisition Time       | 3.2768         |
| 9 Spectrometer Frequency | 500.13         |
| 10 Spectral Width        | 10000.0        |
| 11 Lowest Frequency      | -1914.4        |
| 12 Nucleus               | <sup>1</sup> H |
| 13 Acquired Size         | 32768          |
| 14 Spectral Size         | 65536          |

