



| Parameter | Value |
|--------------------------|---------|
| 1 Solvent | DMSO |
| 2 Temperature | 298.1 |
| 3 Pulse Sequence | zg30 |
| 4 Number of Scans | 8 |
| 5 Receiver Gain | 150 |
| 6 Relaxation Delay | 1.0000 |
| 7 Pulse Width | 8.0000 |
| 8 Acquisition Time | 3.2768 |
| 9 Spectrometer Frequency | 500.13 |
| 10 Spectral Width | 10000.0 |
| 11 Lowest Frequency | -1914.5 |
| 12 Nucleus | 1H |
| 13 Acquired Size | 32768 |
| 14 Spectral Size | 65536 |