


=====

|                                      |   |                                                                             |                   |   |         |
|--------------------------------------|---|-----------------------------------------------------------------------------|-------------------|---|---------|
| Acq. Operator                        | : |                                                                             | Seq. Line         | : | 2       |
| Acq. Instrument                      | : | Instrument 1                                                                | Location          | : | Vial 32 |
| Injection Date                       | : | 3/22/2018 4:41:35 PM                                                        | Inj               | : | 1       |
|                                      |   |                                                                             | Inj Volume        | : | 3.0 µl  |
| Different Inj Volume from Sequence ! |   |                                                                             | Actual Inj Volume | : | 2.0 µl  |
| Acq. Method                          | : | D:                                                                          |                   |   |         |
| Last changed                         | : | 3/24/2017 8:58:09 AM by 1                                                   |                   |   |         |
| Analysis Method                      | : | D:\CHEM32\1\1.M                                                             |                   |   |         |
| Last changed                         | : | 3/22/2018 5:52:27 PM                                                        |                   |   |         |
|                                      |   | (modified after loading)                                                    |                   |   |         |
| Sample Info                          | : | [A = ACN, B = 0.1% H3PO4/H2O] 0% A/B 2 min, 0-95% A/B 5 min, 95% A/B 3 min; |                   |   |         |
|                                      |   | Athena C18 (3 µm, 2.1 mm * 100 mm); Flow: 0.45 ml/min, T = 25               |                   |   |         |

Additional Info : Peak(s) manually integrated



=====  
Area Percent Report  
=====

Sorted By : Signal  
Multiplier : 1.0000  
Dilution : 1.0000  
Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=254,4 Ref=400,100

| Peak # | RetTime [min] | Type | Width [min] | Area [mAU*s] | Height [mAU] | Area %  |
|--------|---------------|------|-------------|--------------|--------------|---------|
| 1      | 6.190         | VV   | 0.0417      | 12.20476     | 4.41540      | 0.2365  |
| 2      | 6.272         | VB   | 0.0374      | 3.88682      | 1.49317      | 0.0753  |
| 3      | 7.334         | BB   | 0.0787      | 9.82920      | 1.71912      | 0.1904  |
| 4      | 7.980         | BB   | 0.0390      | 6.10121      | 2.06058      | 0.1182  |
| 5      | 8.699         | BB   | 0.0475      | 5129.35547   | 1661.48047   | 99.3796 |

Totals : 5161.37746 1671.16873

=====  
\*\*\* End of Report \*\*\*