

| Parameter | Value |
|---------------------------|---------------------|
| 1 Origin | Varian |
| 2 Spectrometer | vnmrs |
| 3 Solvent | dmsd |
| 4 Temperature | 25.0 |
| 5 Pulse Sequence | s2pul |
| 6 Number of Scans | 8 |
| 7 Receiver Gain | 54 |
| 8 Relaxation Delay | 1.0000 |
| 9 Pulse Width | 0.0000 |
| 10 Acquisition Time | 4.0630 |
| 11 Acquisition Date | 2018-06-26T15:31:33 |
| 12 Modification Date | 2018-06-27T06:32:39 |
| 13 Spectrometer Frequency | 399.74 |
| 14 Spectral Width | 8012.8 |
| 15 Lowest Frequency | -1608.0 |
| 16 Nucleus | ¹ H |
| 17 Acquired Size | 32556 |
| 18 Spectral Size | 65536 |

