

Parameter	Value
1 Origin	Varian
2 Spectrometer	vnmrs
3 Solvent	cdcl3
4 Temperature	25.0
5 Pulse Sequence	s2pul
6 Number of Scans	8
7 Receiver Gain	46
8 Relaxation Delay	1.0000
9 Pulse Width	0.0000
10 Acquisition Time	4.0630
11 Acquisition Date	2018-08-08T15:45:41
12 Modification Date	2018-08-09T06:46:40
13 Spectrometer Frequency	399.74
14 Spectral Width	8012.8
15 Lowest Frequency	-1608.0
16 Nucleus	¹ H
17 Acquired Size	32556
18 Spectral Size	65536

