

Parameter	Value
1 Spectrometer	WNMR-I-400MHz
2 Solvent	DMSO
3 Temperature	291.6
4 Relaxation Delay	
5 Pulse Width	
6 Acquisition Time	3.0000
7 Acquisition Date	2019-11-29T11:18:41
8 Modification Date	2019-11-29T11:19:04
9 Spectrometer Frequency	399.95
10 Spectral Width	11013.0
11 Lowest Frequency	-4510.6
12 Nucleus	1H
13 Acquired Size	33039
14 Spectral Size	131072

