

Parameter	
1	Spectrometer WNMN-1-400MHz
2	Solvent DMSO
3	Temperature 291.0
4	Relaxation Delay
5	Pulse Width
6	Acquisition Time 3.0000
7	Acquisition Date 2019-12-03T12:20:01
8	Modification Date 2019-12-03T12:20:24
9	Spectrometer Frequency 399.95
10	Spectral Width 11013.0
11	Lowest Frequency -4510.6
12	Nucleus 1H
13	Acquired Size 33039
14	Spectral Size 131072

