

Parameter	Value
1 Spectrometer	spect
2 Solvent	D2O
3 Temperature	300.0
4 Relaxation Delay	1.0000
5 Pulse Width	13.2236
6 Acquisition Time	3.8011
7 Acquisition Date	2020-03-12T12:35:00
8 Modification Date	2020-03-12T12:35:15
9 Spectrometer Frequency	400.13
10 Spectral Width	8620.7
11 Lowest Frequency	-2409.7
12 Nucleus	<sup>1</sup> H
13 Acquired Size	32768
14 Spectral Size	32768

