

| Parameter                 | Value  |
|---------------------------|--|
| 1 Title                   | 201126.2.fid                                 |
| 2 Origin                  | Bruker BioSpin GmbH                          |
| 3 Owner                   | root   |
| 4 Spectrometer            | spect  |
| 5 Solvent                 | DMSO   |
| 6 Temperature             | 298.1  |
| 7 Pulse Sequence          | zg30   |
| 8 Experiment              | 1D   |
| 9 Probe                   | Z139811_0003 (PA QXI 500SI H/ P-C/ N-D-05 Z) |
| 10 Number of Scans        | 16   |
| 11 Receiver Gain          | 167  |
| 12 Pulse Width            | 10.0000                                      |
| 13 Acquisition Time       | 3.2768                                       |
| 14 Acquisition Date       | 2020-11-26T13:50:53                          |
| 15 Modification Date      | 2020-11-26T13:50:54                          |
| 16 Spectrometer Frequency | 500.13                                       |
| 17 Spectral Width         | 10000.0                                      |
| 18 Lowest Frequency       | -1911.7                                      |
| 19 Nucleus                | <sup>1</sup> H                               |
| 20 Acquired Size          | 32768  |
| 21 Spectral Size          | 65536  |

