



File D:\DATA  
 Injection Date : 09-Mar-15, 15:45  
 Sample Name : 1001264-89-6  
 Acq. Operator : HPLC  
 Spec. Reported : MS Integration  
 Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M  
 Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M  
 Sample Info :  
 Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)  
 Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B  
 Flow Rate: 1.2 ml/min  
 Column: Eclipse XDB-C18,4.6\*150mm,5um  
 Oven Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
7.00	0.07	5351.26	1224.52	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
7.00	0.07	2021.86	464.92	100.00