

File D:\DATA\1025687-58-4 2014060304

Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Sample Info :

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)


Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B

Flow Rate: 1.2 ml/min


Column: Eclipse XDB-C18,4.6*150mm,5um

Oven Temperature: 40 C

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
7.48	0.04	9747.00	3243.36	98.67
8.62	0.04	124.70	52.89	1.26
10.83	0.03	6.29	2.57	0.06

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
7.49	0.04	5927.20	2357.73	98.85
8.62	0.04	68.92	28.97	1.15