


LC/MS Report

File D:\DATA		
Injection Date	: 16 May 14 11:14 am +0800	Tgt Mass(EZX) :
Sample Name	: 1095173-27-5 2014051508	Location : P1-C-01
Acq. Operator	: Product	Inj : 1
Spec. Reported	: MS Integration	Inj Volume : 1 ul
Acq. Method	: D:\METHODS\1-POS-MON.M	
Analysis Method	: D:\METHODS\1-POS-MON.M	
Sample Info : Easy-Access Method: '1-POS-MON'		
Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN		
Gradient: 5%B increase to 95%B within1.2min,		
95%B for 1.3min,back to 5%B within 0.01min.		
Flow Rate: 2 mL/min		
Column: XBridge, 3.5um, 50*4.6mm		
Oven Temperature: 50 C		
Mass Range:110-1000		



Integration Results for ADC1 A, ELSD				
RetTim	Width	Area	Height	Area%
1.27	0.01	1215.33	1293.80	100.00
Integration Results for DAD1 A, Sig=214,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.27	0.01	1215.33	1293.80	100.00
Integration Results for DAD1 B, Sig=254,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.27	0.01	870.80	1155.69	100.00

