


MS ~~File~~D:\DATA
Injection Date : 28 Aug 14 06:45 pm
Sample Name : 1116235-97-2
Acq. Operator : XXXX
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10 mM NH4HCO3) B:ACN
Gradient: 5%B increase to 95%B within1.2min,
95%B for 1.3min,back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: Xbridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range:100-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.59	0.02	1.46	0.92	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.58	0.01	745.64	1052.52	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.58	0.01	511.04	723.55	100.00

MS Report

Ret. Time: 1.58 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.64

