


LC/MS Report

File D:\DATA
 Injection Date : 10 Jul 14 06:45 pm
 Sample Name : 1226056-71-8 Location :
 Acq. Operator : XXXX Inj : 1
 Spec. Reported : MS Integration Inj Volume : 2 ul
 Acq. Method : D:\METHODS\1-POS-MON.M
 Analysis Method : D:\METHODS\1-POS-MON.M
 Sample Info : Easy-Access Method: '1-POS-MON'
 Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
 Gradient: 5%B increase to 95%B within1.2min,
 95%B for 1.3min,back to 5%B within 0.01min.
 Flow Rate: 2 mL/min
 Column: XBridge, 3.5um, 50*4.6mm
 Oven Temperature: 50 C
 Mass Range:110-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.21	0.02	1287.07	1286.70	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.21	0.02	1287.07	1286.70	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.21	0.01	1081.92	1489.28	100.00

LC/MS Report


Ret. Time: 1.19 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.34



Ret. Time: 1.79



Ret. Time: 1.82

