




File D:\DATA
 Injection Date : 28-Jan-16, 20:45
 Sample Name : 1256580-46-7
 Acq. Method : C:\CHEM32\1\METHODS\2-HPLC-15MIN.M
 Analysis Method : C:\CHEM32\1\METHODS\2-HPLC-15MIN.M
 Sample Info : Easy-Access Method: '2-HPLC-15MIN'
 Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B
 Flow Rate: 1.2 ml/min
 Column: Eclipse XDB-C18,4.6*150mm,5um
 Oven Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
8.50	0.06	2793.58	764.04	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
8.50	0.06	1912.01	545.86	100.00