



LC/MS Report

File D:\DATA  
Injection Date : 22-May-14, 11:07:31 Tgt Mass (EZ) :  
Sample Name : 154229-18-2 Location : P1-C-03  
Acq. Operator : HPLC Inj : 0  
Spec. Reported : MS Integration Inj Volume : 10 ul  
Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M  
Analysis Method : C:\CHEM32\1\METHODS\2-HPLC-15MIN.M  
Sample Info :  
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)  
Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B  
Flow Rate: 1.2 ml/min  
Column: Eclipse XDB-C18,4.6\*150mm,5um  
Oven Temperature: 40 C



Integration Results for DAD1 E, Sig=215,4 Ref=off

| RetTim | Width | Area    | Height  | Area%  |
|--------|-------|---------|---------|--------|
| 10.44  | 0.06  | 5269.24 | 1177.45 | 100.00 |

Integration Results for DAD1 B, Sig=254,4 Ref=off

| RetTim | Width | Area    | Height  | Area%  |
|--------|-------|---------|---------|--------|
| 10.44  | 0.08  | 6938.73 | 1543.32 | 100.00 |