



File D:\DATA\
Injection Date : 28 Mar 16 5:11 pm +0800 Tgt Mass (EZX) :
Sample Name : Location :
Acq. Operator : Inj : 1
Spec. Reported : MS Integration Inj Volume : 3 ul
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN'
Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
Mobile phase: H2O(0.05%TFA) (A) / ACN(0.05%TFA) (B)
Elution program: Gradient from 5 to100% of B in 1.3 min at 2ml/min
Temperature: 50°C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 80



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.18	0.14	450.19	39.61	19.23
1.38	0.02	1890.93	1537.79	80.77

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.38	0.02	1006.21	810.47	100.00

Ret. Time: 1.37

<<<< POSITIVE SPECTRA >>>>

