


File D:\DATA\ 162401-32-3 2014030704
 Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
 Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
 Sample Info : *** No target masses specified ***
 Method Info : Mobile Phase: A: water(0..01%TFA) B: ACN(0..01%TFA)
 Gradient: 5%-95% B in 1.5min
 Flow Rate: 2.0ml/min
 Column: Sunfire C18 3.5um,4.6*50mm
 Oven Temperature:50 C
 Mass Range: 110-800 UV: 214nm;254nm



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.92	0.01	2105831.75	2425542.00	100.00


Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.92	0.03	3059.16	1328.77	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.92	0.02	1658.41	1183.45	100.00

Ret. Time: 1.31 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.92



Ret. Time: 2.37

