




File D:\DATA\
 Injection Date : 20 Feb 17 3:19 pm +0800 Tgt Mass(EZX) :
 Sample Name : Location :
 Acq. Operator : Inj : 1
 Spec. Reported : MS Integration Inj Volume : 2 ul
 Acq. Method : D:\METHODS\1-POS-3MIN.M
 Analysis Method : D:\METHODS\1-POS-3MIN.M
 Sample Info : Easy-Access Method: '1-POS-3MIN.M'
 Method Info : Column:SUNFIRE C18 4.6X50mm,3.5um
 Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01% TFA)(B)
 Elution program: Gradient from 5 to 95% of B in1.4min at 2.0ml/min
 Temperature: 50°C
 Detection: UV (214,254 nm) and MS (ESI, POS mode ,113 to 1200 amu)




Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.34	0.02	3014.41	2545.20	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.34	0.02	1382.50	1208.06	100.00

Ret. Time: 1.34 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.50

