





File D:\DATA
Injection Date : 30 Mar 15 3:45
Sample Name : 175131-60-9
Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Sample Info : *** No target masses specified ***
Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN
Gradient: 5%-95% B in 1.5min
Flow Rate: 2.0ml/min
Column: XBridge C18,4.6*50mm,3.5um
Oven Temperature: 40 C
MASS Range:110-800



Integration Results for ELS1 A, Voltage

RetTim	Width	Area	Height	Area%
1.66	0.04	3386.78	1169.68	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.66	0.04	3386.78	1169.68	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.66	0.04	3386.78	1169.68	100.00

Ret. Time: 1.66

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.94

