

File D:\DATA\

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(0..01%TFA) B: ACN(0..01%TFA)

Gradient: 5%-95% B in 1.5min


Flow Rate: 2.0ml/min

Column: Sunfire C18 3.5um,4.6*50mm


Oven Temperature:50 C

Mass Range: 110-800 UV: 214nm;254nm


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.18	0.03	3916377.75	2427535.75	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.17	0.03	2675.70	1315.85	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.17	0.02	2270.69	1659.72	100.00

Ret. Time: 1.16

<<<< POSITIVE SPECTRA >>>>

