



File D:\DATA\
Injection Date : Tgt Mass(EZX) :
Sample Name : Location :
Acq. Operator : Monitor Inj : 1
Spec. Reported : MS Integration Inj Volume : 0.36 u1
Acq. Method : C:\CHEM32\1\METHODS\1-POS-3MIN.M
Analysis Method : C:\CHEM32\1\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN'
Method Info : Mobile Phase: A :Water(0.01%TFA) B:ACN(0.01%TFA)
Gradient: 5%B increase to 95%B within 1.2min,
95%B for 1.3min,back to 5%B within 0.01min
Flow Rate: 2.0ml/min
Column : SunFire C18, 4.6*50mm, 3.5um
Column Temperature: 50 C
Detection:UV(214,4nm) and MS(ESI,POS mode,110 to 1000 amu)


*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.34	0.03	943.28	506.30	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.34	0.03	237.02	127.65	98.82
1.73	0.02	2.83	2.85	1.18

Ret. Time: 1.36


<<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.10



Ret. Time: 2.25



Ret. Time: 2.60

