


LC/MS Report

File D:\DATA
Injection Date : 4 May 14 1:39 pm +0800 Tgt Mass (EZ) :
Sample Name : 212631-79-3 2014050403 Location : P1-C-05
Acq. Operator : Monitor Inj : 1
Spec. Reported : MS Integration Inj Volume : 0.5 ul
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
Gradient: 5%B increase to 95%B within1.2min,
95%B for 1.3min,back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: XBridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range:110-1000



Integration Results for ADC1 A, ELSD				
RetTim	Width	Area	Height	Area%
1.68	0.05	28.54	8.64	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.66	0.01	594.56	772.18	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.66	0.01	219.00	291.44	100.00

Ret. Time: 1.64 <<<< POSITIVE SPECTRA >>>>

