

YX-314045-39-1
YX-314045-39-1

12.63

7.29
7.29
7.28
7.27
7.24
7.24
7.22

3.76
3.29
3.23
3.21
3.19
2.90
2.89
2.81
2.47
2.47
2.46

Parameter	Value
1 Comment	YX-314045-39-1
2 Spectrometer	vnmrs
3 Solvent	dmsd
4 Temperature	25.0
5 Pulse Sequence	s2pul
6 Number of Scans	8
7 Receiver Gain	44
8 Relaxation Delay	1.0000
9 Pulse Width	0.0000
10 Acquisition Time	4.0630
11 Acquisition Date	2017-07-10T16:47:34
12 Modification Date	2017-07-11T07:48:33
13 Spectrometer Frequency	399.74
14 Spectral Width	8012.8
15 Lowest Frequency	-1608.0
16 Nucleus	1H
17 Acquired Size	32556
18 Spectral Size	65536

