


File D:\DATA\38748-32-2 2014030705
 Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
 Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
 Sample Info : *** No target masses specified ***
 Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%-95% B in 1.5min
 Flow Rate: 2.0ml/min
 Column: Sunfire C18 3.5um,4.6*50mm
 Oven Temperature:50 C
 Mass Range: 110-800 UV: 214nm;254nm



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.57	0.02	52857.20	36497.49	100.00


Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.57	0.02	508.76	351.16	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%


Ret. Time: 1.57 <<<< POSITIVE SPECTRA >>>>




Ret. Time: 2.27



Ret. Time: 2.38




Ret. Time: 2.46



LC/MS Report

Ret. Time: 2.54



Ret. Time: 2.76

