


LC/MS Report

File D:\DATA
 Injection Date : 10 Jul 14 06:45 pm
 Sample Name : 391210-10-9
 Acq. Operator : XXXX
 Spec. Reported : MS Integration
 Acq. Method : D:\METHODS\1-POS-MON.M
 Analysis Method : D:\METHODS\1-POS-MON.M
 Sample Info : Easy-Access Method: '1-POS-MON'
 Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
 Gradient: 5%B increase to 95%B within1.2min,
 95%B for 1.3min,back to 5%B within 0.01min.
 Flow Rate: 2 mL/min
 Column: XBridge, 3.5um, 50*4.6mm
 Oven Temperature: 50 C
 Mass Range:110-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.33	0.02	1183.68	1244.58	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.33	0.01	904.71	1187.14	100.00

LC/MS Report

Ret. Time: 1.31 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.76



Ret. Time: 1.89

