

File D:\DATA

Injection Date : 09-Oct-14, 6:45

Sample Name : 41621-49-2

Acq. Operator : HPLC

Spec. Reported : MS Integration

Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Sample Info :

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)


Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B

Flow Rate: 1.2 ml/min


Column: Eclipse XDB-C18,4.6*150mm,5um

Oven Temperature: 40 C

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
7.27	0.04	49.29	20.69	1.14
8.71	0.04	16.02	5.68	0.37
9.73	0.04	4273.93	1502.07	98.49

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
7.08	0.06	14.06	3.45	7.31
7.24	0.04	9.73	3.12	5.06
9.73	0.04	162.53	60.29	84.50
10.41	0.04	6.04	2.37	3.14