


260632 H1-NMR CDC13 303K AV-500




8.0463  
7.7500  
7.4660  
7.4484  
7.3796  
7.3771  
7.3620  
7.3595  
7.2579

0.0000

NAME km1  
EXPNO 170  
PROCNO 1  
Date\_ 20121228  
Time 16.27  
INSTRUM av500  
PROBHD 5 mm QNP 1H/15  
PULPROG zg30  
TD 32768  
SOLVENT CDC13  
NS 1  
DS 0  
SWH 10000.000 Hz  
FIDRES 0.305176 Hz  
AQ 1.6385000 sec  
RG 32  
DW 50.000 usec  
DE 6.50 usec  
TE 303.0 K  
D1 1.00000000 sec  
TD0 1

===== CHANNEL f1 =====  
NUC1 1H  
P1 14.80 usec  
PL1 -1.00 dB  
PL1W 14.32511806 W  
SFO1 500.1340010 MHz  
SI 32768  
SF 500.1300149 MHz  
WDW EM  
SSB 0  
LB 1.00 Hz  
GB 0  
PC 1.00



1.000  
1.023  
1.058  
1.046