



File D:\DATA

Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M
 Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M
 Sample Info : *** No target masses specified ***
 Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B
 Flow Rate: 1.2 ml/min
 Column: Eclipse XDB-C18,4.6*150mm,5um
 Oven Temperature: 40 C


ADC1 A, ELSD



DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
11.70	0.04	90.65	30.78	2.83
12.13	0.20	791.64	47.40	24.70
12.31	0.04	152.58	54.53	4.76
12.36	0.02	152.19	90.83	4.75
12.38	0.02	90.07	68.76	2.81
12.42	0.03	205.83	89.58	6.42
12.47	0.05	259.57	73.46	8.10
12.58	0.02	34.75	23.67	1.08
12.74	0.08	342.49	56.55	10.69
13.01	0.02	20.48	17.96	0.64
13.06	0.05	223.11	52.46	6.96
13.36	0.02	53.71	36.76	1.68
13.43	0.05	94.46	24.34	2.95
13.61	0.02	27.45	24.98	0.86
13.75	0.10	349.47	42.93	10.90
14.12	0.05	66.29	20.78	2.07
14.20	0.04	100.78	29.40	3.14
14.42	0.00	-0.22	6.82	-0.01
14.84	0.05	108.07	29.26	3.37
14.99	0.04	41.72	15.92	1.30

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
9.34	0.05	10147.98	2527.69	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
9.34	0.05	689.85	226.08	100.00