





LC/MS Report

File D:\DATA
Injection Date : 11 Jun 14 06:45 pm
Sample Name : 501-36-0 2014061101
Acq. Operator : xxxxxx
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10mMNH₄HCO₃) B: ACN
Gradient: 5%B increase to 95%B within 1.2min,
95%B for 1.3min, back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: XBridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range: 110-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.21	0.01	2155.05	2521.87	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.21	0.01	2155.05	2521.87	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.21	0.01	412.87	514.93	100.00

Ret. Time: 1.20 <<<< POSITIVE SPECTRA >>>>

