




File D:\DATA\
 Injection Date : 3 Mar 17 1:10 pm +0800 Tgt Mass(EZX) :
 Sample Name : Location :
 Acq. Operator : Inj : 1
 Spec. Reported : MS Integration Inj Volume : 1 ul
 Acq. Method : D:\METHODS\1-POS-3MIN.M
 Analysis Method : D:\METHODS\1-POS-3MIN.M
 Sample Info : Easy-Access Method: '1-POS-3MIN.M'
 Method Info : Column:SUNFIRE C18 4.6X50mm,3.5um
 Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01% TFA)(B)
 Elution program: Gradient from 5 to 95% of B in1.4min at 2.0ml/min
 Temperature: 50°C
 Detection: UV (214, 254 nm) and MS (ESI, POS mode ,113 to 1200 amu)




Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.44	0.02	856.68	815.38	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.44	0.02	1149.51	1092.42	100.00


Ret. Time: 1.46 <<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.74




Ret. Time: 2.33




Ret. Time: 2.38




Ret. Time: 2.48




Ret. Time: 2.71



Ret. Time: 2.75



Ret. Time: 2.82



Ret. Time: 2.90

