




```

File D:\DATA\
Injection Date   : 12 Jun 16 12:19 pm +0800      Tgt Mass(EZX) :
Sample Name     :                               Location    :
Acq. Operator  :                               Inj         : 1
Spec. Reported  : MS Integration                Inj Volume  : 1 ul
Acq. Method    : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info    : Easy-Access Method: '1-POS-3MIN.M'
Method Info    : Column:SUNFIRE C18 4.6X50mm,3.5um
                  Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01% TFA)(B)
                  Elution program: Gradient from 5 to 95% of B in1.4min at 2.0ml/min
                  Temperature: 45°C
                  Detection: UV (214 , 4 nm) and MS (ESI, POS mode ,113 to 1200 amu)
    
```



Integration Results for DAD1 A, Sig=214,4 Ref=off

| RetTim | Width | Area | Height | Area% |
|--------|-------|--------|--------|--------|
| 1.50 | 0.01 | 442.54 | 564.67 | 100.00 |

Integration Results for DAD1 B, Sig=254,4 Ref=off

| RetTim | Width | Area | Height | Area% |
|--------|-------|--------|--------|--------|
| 1.50 | 0.01 | 365.61 | 467.01 | 100.00 |

Ret. Time: 1.50

<<<< POSITIVE SPECTRA >>>>

