




```

File D:\DATA\
Injection Date   : 29 Mar 16   5:02 pm +0800       Tgt Mass(EZX) :
Sample Name     :                               Location   :
Acq. Operator   :                               Inj       : 1
Spec. Reported  : MS Integration                 Inj Volume  : 1 ul
Acq. Method     : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info    : Easy-Access Method: '1-POS-3MIN.M'
Method Info    : Column:SUNFIRE C18 4.6X50mm,3.5um
                  Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01% TFA)(B)
                  Elution program: Gradient from 5 to 95% of B in1.4min at 2.0ml/min
                  Temperature: 45°C
                  Detection: UV (214 , 4 nm) and MS (ESI, POS mode ,113 to 1200 amu)
    
```



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.77	0.01	730.91	876.28	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.77	0.01	1858.47	2234.73	100.00

Ret. Time: 1.77

<<<< POSITIVE SPECTRA >>>>

