


LC/MS Report

File D:\DATA
Injection Date : 16 Jun 14 06:45 pm
Sample Name : 6146-52-7 2014061611 Location :
Acq. Operator : XXXX Inj : 1
Spec. Reported : MS Integration Inj Volume : 1 ul
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
Gradient: 5%B increase to 95%B within1.2min,
95%B for 1.3min,back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: XBridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range:110-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.41	0.02	1227.50	1254.66	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.41	0.01	1627.59	1746.94	100.00

LC/MS Report

Ret. Time: 1.39 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.43

