

YX-659730-32-2
YX-659730-32-2

Parameter	Value
1 Spectrometer	vnmrs
2 Solvent	dmsd
3 Temperature	25.0
4 Pulse Sequence	s2pul
5 Number of Scans	8
6 Receiver Gain	54
7 Relaxation Delay	1.0000
8 Pulse Width	0.0000
9 Acquisition Time	4.0630
10 Acquisition Date	2017-07-17T16:21:04
11 Spectrometer Frequency	399.74
12 Spectral Width	8012.8
13 Lowest Frequency	-1608.0
14 Nucleus	1H
15 Acquired Size	32556
16 Spectral Size	65536

