



File D:\DATA\
 Injection Date : 6 Mar 17 3:14 pm +0800 Tgt Mass(EZX) :
 Sample Name : Location :
 Acq. Operator : Inj : 1
 Spec. Reported : MS Integration Inj Volume : 1 ul
 Acq. Method : D:\METHODS\1-POS-3MIN.M
 Analysis Method : D:\METHODS\1-POS-3MIN.M
 Sample Info : Easy-Access Method: '1-POS-3MIN.M'
 Method Info : Column:SUNFIRE C18 4.6X50mm,3.5um
 Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01% TFA)(B)
 Elution program: Gradient from 5 to 95% of B in1.4min at 2.0ml/min
 Temperature: 50°C
 Detection: UV (214,254 nm) and MS (ESI, POS mode ,113 to 1200 amu)


*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70




Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.35	0.02	14.61	12.67	100.00


Ret. Time: 0.38 <<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.37




Ret. Time: 1.45




Ret. Time: 1.49




Ret. Time: 1.59




Ret. Time: 1.65




Ret. Time: 1.79




Ret. Time: 1.92




Ret. Time: 2.08




Ret. Time: 2.27




Ret. Time: 2.34




Ret. Time: 2.46




Ret. Time: 2.65




Ret. Time: 2.71




Ret. Time: 2.78



Ret. Time: 2.83



Ret. Time: 2.88



Ret. Time: 2.92

