


MS Report

Injection Date : 11 Sep 14 06:45  
 Sample Name : 742112-33-0  
 Acq. Operator : xxxx  
 Spec. Reported : MS Integration  
 Acq. Method : D:\METHODS\1-POS-MON.M  
 Analysis Method : D:\METHODS\1-POS-MON.M  
 Sample Info : Easy-Access Method: '1-POS-MON'  
 Method Info : Mobile Phase: A: water(10 mM NH4HCO3) B:ACN  
 Gradient: 5%B increase to 95%B within1.2min,  
 95%B for 1.3min,back to 5%B within 0.01min.  
 Flow Rate: 2 mL/min  
 Column: Xbridge, 3.5um, 50\*4.6mm  
 Oven Temperature: 50 C  
 Mass Range: 100-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.85	0.03	1.27	0.69	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.83	0.01	535.54	570.75	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.83	0.01	836.55	903.88	100.00

MS Report

Ret. Time: 1.83

<<<< POSITIVE SPECTRA >>>>

