


LC/MS Report

File D:\DATA
Injection Date : 13 Oct 14 6:45
Sample Name : 7689-03-4
Acq. Operator : XXXX
Spec. Reported : MS Integration
Acq. Method : C:\Chem32\1\METHODS\1-POS-2MIN.M
Analysis Method : C:\Chem32\1\METHODS\1-POS-2MIN.M
Sample Info : Easy-Access Method: '1-POS-2MIN'
Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%B increase to 100%B within 0.8min,100%B for 1.1min.
Flow Rate:1.5mL/min
Column: Chromolith Fast Gradien RP-18e,50-3mm.
Column Temperature: 40C
Mass Range:100-1100



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
0.85	0.01	371.49	415.52	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.85	0.01	371.49	415.52	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
0.85	0.01	282.94	314.49	100.00

Ret. Time: 0.86 <<<< POSITIVE SPECTRA >>>>

