



File D:\DATA
Injection Date : 24-Dec-14, 15:45
Sample Name : 81129-83-1
Acq. Operator : HPLC
Spec. Reported : MS Integration
Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M
Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M
Sample Info :
Method Info : Mobile Phase: A: water (0.01%TFA) B: ACN (0.01%TFA)
Gradient: 0min 5%B, 3min 5%B, 10min 95%B, 15min 95%B
Flow Rate: 1.2 ml/min
Column: Eclipse XDB-C18, 4.6*150mm, 5um
Oven Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
6.34	0.02	3.20	2.45	0.17
6.36	0.02	3.65	2.42	0.20
6.82	0.03	1805.81	882.61	98.31
7.97	0.04	17.40	4.99	0.95
10.80	0.03	6.84	2.58	0.37

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
6.82	0.03	561.66	275.81	98.79
7.98	0.03	6.87	2.40	1.21