


MS Report

File D:\DATA
 Injection Date : 1 Sep 14 06:45
 Sample Name : dashatini
 Acq. Operator : XXXX
 Spec. Reported : MS Integration
 Acq. Method : D:\METHODS\1-POS-MON.M
 Analysis Method : D:\METHODS\1-POS-MON.M
 Sample Info : Easy-Access Method: '1-POS-MON'
 Method Info : Mobile Phase: A: water(10 mM NH4HCO3) B:ACN
 Gradient: 5%B increase to 95%B within1.2min,
 95%B for 1.3min,back to 5%B within 0.01min.
 Flow Rate: 2 mL/min
 Column: Xbridge, 3.5um, 50*4.6mm
 Oven Temperature: 50 C
 Mass Range:100-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.47	0.03	12.72	7.48	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.33	0.01	24.89	43.11	1.29
1.45	0.01	1899.93	2325.77	98.71


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.45	0.01	1472.84	2357.46	100.00

MS Report

Ret. Time: 1.33

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.45

