

File D:\DATA

Injection Date : 29 Feb 16 12:45

Sample Name : 865854--05-3

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : Easy-Access Method: '1-POS-3min'

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.2min


Flow Rate: 2.2 ml/min

Column: Poroshell 120 EC-C18, 4.6*30mm, 2.7um


Oven Temperature: 50 C

Mass Range: 110-1500


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.76	0.02	8.09	5.36	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.76	0.02	1822.49	1273.52	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.76	0.02	163.83	118.02	100.00


LC/MS Report

Ret. Time: 1.76

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.15



Ret. Time: 2.48

