

YX-86639-52-3
YX-86639-52-3

10.26

8.00
7.98

7.39
7.36
7.21

6.44

5.38
5.24

3.28
3.06
3.04
3.03

2.48
2.47
2.46
2.46

1.85
1.83
1.81
1.79

1.29
1.27
1.25
0.86
0.84
0.82

Parameter	Value
1 Spectrometer	vnmrs
2 Solvent	dmsd
3 Temperature	25.0
4 Pulse Sequence	s2pul
5 Number of Scans	8
6 Receiver Gain	52
7 Relaxation Delay	1.0000
8 Pulse Width	0.0000
9 Acquisition Time	4.0630
10 Acquisition Date	2017-07-17T16:18:13
11 Spectrometer Frequency	399.74
12 Spectral Width	8012.8
13 Lowest Frequency	-1608.0
14 Nucleus	1H
15 Acquired Size	32556
16 Spectral Size	65536

