


LC/MS Report

File D:\DATA
Injection Date : 26 Jun 14 06:45 pm
Sample Name : 868539-98-4
Acq. Operator : XXXX
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
Gradient: 5%B increase to 95%B within1.2min,
95%B for 1.3min,back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: XBridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range:110-1000



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.45	0.05	269.84	85.90	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.43	0.01	761.91	823.94	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.43	0.01	17.42	20.67	100.00

LC/MS Report


Ret. Time: 0.82 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.33



Ret. Time: 1.35



Ret. Time: 1.42

