


LC/MS Report

File D:\DATA
Injection Date : 10 Jul 14 06:45 pm
Sample Name : 956697-53-3
Acq. Operator : XXXX
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-MON.M
Analysis Method : D:\METHODS\1-POS-MON.M
Sample Info : Easy-Access Method: '1-POS-MON'
Method Info : Mobile Phase: A: water(10mMNH4HCO3) B: ACN
Gradient: 5%B increase to 95%B within1.2min,
95%B for 1.3min,back to 5%B within 0.01min.
Flow Rate: 2 mL/min
Column: XBridge, 3.5um, 50*4.6mm
Oven Temperature: 50 C
Mass Range:110-1000




Integration Results for ADC1 A, ELSD				
RetTim	Width	Area	Height	Area%
1.48	0.01	24.01	30.05	1.97
1.64	0.02	1193.56	1227.72	98.03

Integration Results for DAD1 A, Sig=214,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.48	0.01	24.01	30.05	1.97
1.64	0.02	1193.56	1227.72	98.03

Integration Results for DAD1 B, Sig=254,4 Ref=off				
RetTim	Width	Area	Height	Area%
1.64	0.01	855.71	1061.92	100.00

LC/MS Report

Ret. Time: 1.46 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.62

